



## System Requirements

### Hardware

Compass requires a minimum performance graphics card be used in the application PC on which Compass is installed. The following table lists the supported graphics cards. Please contact the Support team for an up-to-date list and information on video card driver versions.

Supported Graphics Cards:

- ATI Radeon 9800 Pro
- ATI Radeon x600
- Quadro FX 1300
- Quadro 980 xgl
- Quadro FX 1500

OS/Language: Windows platform - XP Pro English or greater.

CPU Architecture: Single Core, Dual Core, Quad Core.

### Compliance

This product has been developed in compliance with ANSI/AAMI/ISO 14971:2003, Medical Devices – Application of risk management to medical devices.

### Software

- Requires Windows XP SP3 / Vista / Windows 7
- Microsoft .Net Compact Framework 3.5 is required (installed by installer if not preset)
- Compass uses Microsoft SQL Server Express 2008 (optionally installed by installer)
- SQL Server Instance Named localhost\SQLEXPRESS (optionally created by installer)
- TCP/IP protocol enabled on SQLEXPRESS for port 1433 (optionally configured by installer)
- Shared Host Service (installed as a Windows service with automatic startup by installer)

Compass Launch API Service (installed as windows application, launched at system logon via registry settings)